

GIMA VITAL PRO MONITOR

The Vital Pro Monitor is a multi-parameter patient monitor that measures the following physiological parameters: **non-invasive blood pressure measurement (NIBP)**, **oxygen saturation (SpO₂)**, electrocardiography (ECG) and temperature (TEMP).

Properties:

3.5" LCD display with **large numbers**
 Record display in waveform, list and trend graph
 USB port for data transfer and **software upgrade**

Technical specifications:

Blood pressure Measurement Method: Oscillometric
 Messmodi: Manuel, Auto, STAT
 Messbereich: SYS mmHg: (Adult 40-275) (Pediatric 40-200) (Neonate 40-135)
 MAP mmHg: (Adult 20-230) (Pediatric 20-165) (Neonate 20-110)
 DIA mmHg: (Adult 10-210) (Pediatric 10-150) (Neonate 10-95)

SpO₂ Measuring range: 35% ~ 100% ± 3% (between 70 to 100%)

Low perfusion power: < 0.4%

Pulse 30 ~ 240BPM ± 2 BPM / 2%

ELECTROCARDIOGRAM Display heartbeat: 15~130 bpm ± 1% or ±2 bpm

Alarm if no heartbeat within 10s

Sensitivity: x1/4, x1/2, x1, x2

Speed: 25 mm/s

Standard accessories:

Adult SpO₂ Fingertip, Blood Pressure Cuff for Adults, Temperature Probe, ECG Cable + 10 ECG Electrodes, Lithium Ion Battery (11.1V,4400mAh), Power Cord (100-240V, 50/60Hz), Dust Jacket, User manual (GB, FR, IT, ES, D (soon))

Optional accessories:

Equipment trolley



Item no.: 35124



Item no.: 35127

MEDICARE AG

Hauptstrasse 51, 5024 Küttigen, Phone 044 482 482 6 / Fax 044 482 74 88

info@medicareag.ch / www.medicareag.ch